

# IL26 Rabbit mAb

货号: **AYM30781**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	
实测分子量	40kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse heart
细胞定位	cytosol,extracellular region,extracellular space
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human IL26.
序列	KHKQSSFTKSCYPRGTLTSAVDALYIKAAWLKATIPEDRIKNIRLLKKTKKQFMKNCQFQEQLLSFFMEDVFGQLQLQ

## 靶点信息

研究背景	This gene was identified by its overexpression specifically in herpesvirus samimiri-transformed T cells. The encoded protein is a member of the IL10 family of cytokines. It is a secreted protein and may function as a homodimer. This protein is thought to contribute to the transformed phenotype of T cells after infection by herpesvirus samimiri. [provided by RefSeq, Jul 2008]
基因ID	55801
基因名	IL26

Swiss	Q9NPH9
别名	AK155;IL-26;IL26

## 产品验证



Western blot analysis of IL26 expressed in Mouse heart using IL26 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)