

— ABLYBIO, Help Your Research



Cathepsin L/V/K/H Rabbit mAb

货号: **AYM30852**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	
实测分子量	38kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	
细胞定位	
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Cathepsin L/V/K/H.
------	--

靶点信息

研究背景	
------	--

基因ID	
基因名	
Swiss	P07711O60911P43235P09668 (https://www.uniprot.org/uniprotkb/P07711O60911P43235P09668/entry)
别名	Cathepsin L/V/K/H,Cathepsin L/V/K/H Rabbit mAb

产品验证



Western blot analysis of Cathepsin L/V/K/H expressed in human liver cancer, mouse liver, rat muscle using Cathepsin L/V/K/H Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)