

Menin Rabbit mAb

货号: **AYM30862**

产品信息

反应	Human
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50
理论分子量	63kDa/67kDa/68kDa
实测分子量	68kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,HepG2,K-562
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Menin.
序列	QKLLWLLYDLGHLERYPMALGNLADLEELEPTGPRPDPLTYLHKGIASAKTYRDEHIYPYMYLAGYHCRNRNVREALQAW ADTATVIQDYNYCREDEEIIYKEFFEVDVNDVIPNLLKEAASLLEAGEERPGEQSQGTQSQGSALQDPECFAHLLRFYDGICK WEEGSPTPVLHVGWATFLVQSLGRFEGQVRQKVRIVSREAEAAEAEPEPWGEEAREGRRRGPRRESKPEEPPPPKPPALDK GLGTGQGAVSGPPRPPGTVAGTARGPEGGSTAQVPAPAASPPPEGPVLTQSEKMKGMKELLVATKINSSAIKLQLTAQ SQVQMKKQKVSTPSDYTL SFLKRQRKGL

靶点信息

研究背景	This gene encodes menin, a putative tumor suppressor associated with a syndrome known as multiple endocrine neoplasia type 1. In vitro studies have shown menin is localized to the nucleus, possesses two functional nuclear localization signals, and inhibits transcriptional activation by JunD, however, the function of this protein is not known. Two messages have been detected on northern blots but the larger message has not been characterized. Alternative splicing results in multiple transcript variants.
基因ID	4221
基因名	MEN1
Swiss	O00255
别名	MEN1;MEAI;SCG2;menin

产品验证



Western blot analysis of Menin expressed in HeLa, HepG2, K-562 using Menin Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn