

SULT2A1 Rabbit mAb

货号: **AYM30908**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	33kDa
实测分子量	34kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse skeletal muscle,Mouse liver
细胞定位	Cytoplasm
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human SULT2A1.
序列	MSDDFLWFEGIAFPTMGFRSETLRKVRDEFVIRDEDVILTYPKSGTNWLAIEILCLMHSKGDAKWIQSVPIWERSPPWVESEI GYTALSETESPRLFSSHLPIQLFPKSFSSKAKVIYLMRNPRDVLVSGYFFWKNMKFIKKPKS

靶点信息

研究背景	This gene encodes a member of the sulfotransferase family. Sulfotransferases aid in the metabolism of drugs and endogenous compounds by converting these substances into more hydrophilic water-soluble sulfate conjugates that can be easily excreted. This protein catalyzes the sulfation of steroids and bile acids in the liver and adrenal glands, and may have a role in the inherited adrenal androgen excess in women with polycystic ovary syndrome.
------	--

基因ID	6822
基因名	SULT2A1
Swiss	Q06520
别名	SULT2A1;DHEA-ST;DHEAS;HST;ST2;ST2A1;ST2A3;STD;hSTa

产品验证



Western blot analysis of SULT2A1 expressed in Mouse skeletal muscle, Mouse liver using SULT2A1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn