

# Skp1 Rabbit mAb

货号: **AYM31047**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IF/ICC:</b> 1:50 - 1:200 <b>FC:</b> 1:20 - 1:50
理论分子量	18kDa
实测分子量	19kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	A-549,A-431,HeLa,LO2,SW480,K-562,Mouse kidney
细胞定位	centrosome,cytoplasm,cytosol,nucleoplasm,nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Skp1.
序列	MPSIKLQSSDGEIFEVDVEIAKQSVTIKTMLEDLGMDDDEGDDDPVPLPNVNAAILKKVIQWCTHHKDDPPPPEDDENKEK RTDDIPVWDQEFLKVDQGLFELILAANYLDIKGLLDVTCKTVANMIKGGKTPEEIRKTFNIKNDFTEEEEAQVRKENQWCEE K

## 靶点信息

研究背景	This gene encodes a component of SCF complexes, which are composed of this protein, cullin 1, a ring-box protein, and one member of the F-box family of proteins. This protein binds directly to the F-box motif found in F-box proteins. SCF complexes are involved in the regulated ubiquitination of specific protein substrates, which targets them for degradation by the proteasome. Specific F-box proteins recognize different target protein(s), and many specific SCF substrates have been identified including regulators of cell cycle progression and development. Studies have also characterized the protein as an RNA polymerase II elongation factor. Alternative splicing of this gene results in two transcript variants. A related pseudogene has been identified on chromosome 7.
基因ID	6500
基因名	SKP1
Swiss	P63208
别名	SKP1;EMC19;OCP-II;OCP2;SKP1A;TCEB1L;p19A

## 产品验证



Western blot analysis of Skp1 expressed in Mouse kidney, A-549, A-431 using Skp1 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)