

Hsp90 beta Rabbit mAb

货号: **AYM31225**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	83kDa
实测分子量	85kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	K-562,DU145,SW620,COS-7,NIH/3T3,HeLa
细胞定位	Cytoplasm,Melanosome
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Hsp90 beta.
序列	MPEEVHHGEEEVETFAFQAEIAQLMSLIINTFYSNKEIFLRELISNASDALDKIRYESLTDP SKLDSGKELKIDIIPNPQERTLTLVDTGIGMTKADLNNLGTIAKSGTKAFMEALQAGADISMIGQFGVGFYSAYLVAEKVVVITKHNDDEQYAWESSAGGSFTVRADHGEPGRGTVLHLKEDQTEYLEERRVKEVVKKHSQFIGYPITLYLEKEREKEISDDEAEKEEKEEEDKDDEEKPKIEDVGSDEEDDSGKDKKKKTKKIKEYIDQEELNKTPIWTRNPDDITQEE


靶点信息

研究背景	This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes.
基因ID	3326
基因名	HSP90AB1
Swiss	P08238
别名	D6S182;HSP84;HSP90B;HSPC2;HSPCB;Hsp90 beta;HSP90AB1

产品验证



Immunohistochemical analysis of paraffin-embedded human testis, using Hsp90 beta Antibody.



Western blot analysis of Hsp90 beta expressed in K-562,DU145,SW620 using Hsp90 beta Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn