

— ABLYBIO, Help Your Research



Cdc34 Rabbit mAb

货号: **AYM31227**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IF/ICC IP FC
推荐浓度	WB: 1:500 - 1:2000 IF/ICC: 1:50 - 1:200 IP: 1:20 - 1:50 FC: 1:20 - 1:50
理论分子量	26kDa
实测分子量	35kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Jurkat,C6,SH-SY5Y,Mouse testis
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein corresponding to Human Cdc34.
------	--

靶点信息

研究背景	The protein encoded by this gene is a member of the ubiquitin-conjugating enzyme family. Ubiquitin-conjugating enzyme catalyzes the covalent attachment of ubiquitin to other proteins. This protein is a part of the large multiprotein complex, which is required for ubiquitin-mediated degradation of cell cycle G1 regulators, and for the initiation of DNA replication.
基因ID	997
基因名	CDC34
Swiss	P49427
别名	Cdc34,Cdc34 Rabbit mAb,(E3-independent) E2 ubiquitin-conjugating enzyme R1,E2 ubiquitin-conjugating enzyme R1,Ubiquitin-conjugating enzyme E2-32 kDa complementing,Ubiquitin-conjugating enzyme E2-C DC34,Ubiquitin-protein ligase R1,UBCH3,UBE2R1

产品验证



Western blot analysis of Cdc34 expressed in Mouse testis, HeLa, Jurkat using Cdc34 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)