

— ABLYBIO, Help Your Research



Hsp90 beta Rabbit mAb

货号: **AYM31336**

产品信息

| | |
|-------|---|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Monoclonal |
| 预测反应 | |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 |
| 理论分子量 | 83kDa |
| 实测分子量 | 90kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa cells |
| 细胞定位 | Cytoplasm,Melanosome |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein corresponding to Human Hsp90 beta . |
|------|--|


靶点信息

| | |
|-------|--|
| 研究背景 | This gene encodes a member of the heat shock protein 90 family; these proteins are involved in signal transduction, protein folding and degradation and morphological evolution. This gene encodes the constitutive form of the cytosolic 90 kDa heat-shock protein and is thought to play a role in gastric apoptosis and inflammation. Alternative splicing results in multiple transcript variants. Pseudogenes have been identified on multiple chromosomes. |
| 基因ID | 3326 |
| 基因名 | HSP90AB1 |
| Swiss | P08238 |
| 别名 | Hsp90 beta,Hsp90 beta Rabbit mAb,HSP90AB1,Heat shock 84 kDa,Heat shock protein family C member 3,HSP90B,HSPC2,HSPC3,HSPCB |

产品验证



Western blot analysis of Hsp90 beta expressed in HeLa cells using Hsp90 beta Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.



Immunohistochemical analysis of paraffin-embedded human kidney, using Hsp90 beta antibody .

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)