

Exportin 5 Rabbit mAb

货号: **AYM31375**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC FC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 FC: 1:20 - 1:50
理论分子量	136kDa
实测分子量	136kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF7,U-251MG,NIH/3T3,HeLa,PC-12,Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Exportin 5.
序列	MRLSQKWQVINQRSLLCGEDEAADENPESQEMLEEQLVRLMLTREVMDLITVCCVSKKGADHSSAPPADGDDEEMMATE VTPSAMAELTDLGKCLMKHEDVCTALLITAFNSLAWKDTLSCQRRTTSQLCWPLLKQVLSGTLADAVTWLFTSVLKGLQM HGQHDGCMASLVHLAFQIYEALRPRYLEIRAVMEQIPEIQKDSLQDQFCKLLNPSLQKVADKRRKQDFKRLIAGCIGKPLGE QFRKEVHIKNLPSLFKKTTPMLETEVLNDGGGLATIFEP

靶点信息

研究背景	This gene encodes a member of the karyopherin family that is required for the transport of small RNAs and double-stranded RNA-binding proteins from the nucleus to the cytoplasm. The encoded protein translocates cargo through the nuclear pore complex in a RanGTP-dependent process.
------	--

基因ID	57510
基因名	XPO5
Swiss	Q9HAV4
别名	XPO5;exp5

产品验证



Western blot analysis of Exportin 5 expressed in Mouse brain,Rat brain,MCF7 using Exportin 5 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30u g per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn