

# Hsp90 Rabbit mAb

货号: **AYM31393**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Monoclonal
预测反应	
应用	WB IHC IF/ICC IP FC
推荐浓度	<b>WB:</b> 1:500 - 1:2000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:20 - 1:50 <b>FC:</b> 1:20 - 1:50
理论分子量	90kDa
实测分子量	90kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.75% BSA,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,Raji,Mouse liver,Rat liver,Rat brain
细胞定位	Cell membrane,Cytoplasm,Melanosome
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein corresponding to Human Hsp90 .
序列	LKEFEGKTLVSVTKEGLELPEDEEEKQEEKTKFENLCKIMKDILEKKVEKVVVSNRLVTSPPCCIVTSTYGWTANMERIMK AQLRDNSTMGYMAAKKHLEINPDHSIIETLRQKAEADKNDKSVKDLVILLYETALLSSGFSLEDQPQTHANRIYRMIKLGLGI DEDDPTADDTSAAVTEEMPPLEGDDDDTSRMEVD

## 靶点信息


研究背景	The protein encoded by this gene is an inducible molecular chaperone that functions as a homodimer. The encoded protein aids in the proper folding of specific target proteins by use of an ATPase activity that is modulated by co-chaperones. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2012]
------	---

基因ID	3320,3326
基因名	HSP90AA1,HSP90AB1
Swiss	P07900/P08238
别名	EL52; HEL-S-65p; HSP86; HSP89A; HSP90A; HSP90N; HSPC1; HSPCA; HSPCAL1; HSPCAL4; HSPN; Hsp103; Hsp89; Hsp90; LAP-2; LAP2

## 产品验证



Immunohistochemical analysis of paraffin-embedded Rat spleen, using the Antibody at 1:300 dilution.



Western blot analysis of Hsp90 expressed in Mouse liver, Rat liver, HeLa using Hsp90 Rabbit mAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)