

Lsh/HELLS Rabbit pAb

货号: **AYP10299**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IP
推荐浓度	WB: 1:500 - 1:2000 IP: 1:50 - 1:200
理论分子量	14kDa/36-39kDa/81-97kDa
实测分子量	120kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	MCF-7,A549,293T,HepG2,Mouse thymus
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 719-838 of human Lsh/HELLS (NP_060533.2).
序列	VYRLVTANTIDQKIVERAAAKRKLEKLIHKNHFKGGQSGLNLSKNFLDPKELMELLKSRDYEREIKGSREKVISDKDLELLLD RSDLIDQMNASGPIKEKMGIFKILENSEDSSPECLF

靶点信息

研究背景	This gene encodes a lymphoid-specific helicase. Other helicases function in processes involving DNA strand separation, including replication, repair, recombination, and transcription. This protein is thought to be involved with cellular proliferation and may play a role in leukemogenesis. Alternatively spliced transcript variants encoding different isoforms have been identified.
基因ID	3070

基因名	HELLS
Swiss	Q9NRZ9
别名	HELLS;ICF4;LSH;Nbla10143;PASG;SMARCA6

产品验证



Western blot analysis of Lsh/HELLS expressed in Mouse thymus, MCF-7, A549 using Lsh/HELLS Rabbit p Ab at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 12 0s.

实验步骤

访问官网浏览详情: www.ablybio.cn