

CSNK1E Rabbit pAb

货号: **AYP10313**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens
应用	WB IHC IF/ICC IP
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200 IP: 1:50 - 1:100
理论分子量	47kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 247-416 of human CSN K1E (NP_689407.1).
序列	EFSTYLNFCRSLRFDDKPDYSYLRQLFRNLFHRQGFSDYVFDWNMLKFGAARNPEDVDREHEREEREMGQLRGSAT RALPPGPPTGATANRLRSAAEPVASTPASRIQPAGNTSPRAISRVDREKVSMLRHRGAPANVSSSDLTGRQEVSRIPASQ TSVPFDHLGK

靶点信息

研究背景	The protein encoded by this gene is a serine/threonine protein kinase and a member of the casein kinase I protein family, whose members have been implicated in the control of cytoplasmic and nuclear processes, including DNA replication and repair. The encoded protein is found in the cytoplasm as a monomer and can phosphorylate a variety of proteins, including itself. This protein has been shown to phosphorylate period, a circadian rhythm protein. Two transcript variants encoding the same protein have been found for this gene.
基因ID	1454
基因名	
Swiss	P49674
别名	CSNK1E;CKIepsilon;HCKIE

产品验证



Western blot analysis of CSNK1E expressed in HeLa,293T using CSNK1E Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn