

— ABLYBIO, Help Your Research



# MST1/STK4 Rabbit pAb

货号: **AYP10388**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IF:</b> Homo sapiens <b>WB:</b> Mus musculus
应用	WB IHC IF/ICC IP
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:50 - 1:200
理论分子量	52kDa/55kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,U-87MG,HepG2,HeLa,MCF-7,Mouse brain,Mouse lung,Mouse spleen,Rat liver
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 303-487 of human MST 1/STK4 (NP_006273.1).
------	--

## 靶点信息

研究背景	The protein encoded by this gene is a cytoplasmic kinase that is structurally similar to the yeast Ste20p kinase, which acts upstream of the stress-induced mitogen-activated protein kinase cascade. The encoded protein can phosphorylate myelin basic protein and undergoes autophosphorylation. A caspase-cleaved fragment of the encoded protein has been shown to be capable of phosphorylating histone H2B. The particular phosphorylation catalyzed by this protein has been correlated with apoptosis, and it's possible that this protein induces the chromatin condensation observed in this process.
基因ID	6789
基因名	STK4
Swiss	Q13043 ( <a href="https://www.uniprot.org/uniprotkb/Q13043/entry">https://www.uniprot.org/uniprotkb/Q13043/entry</a> )
别名	KRS2,MST1,YSK3,STK4,MST1/STK4 Rabbit pAb,Mammalian STE20-like protein kinase 1,STE20-like kinase MST1,Serine/threonine-protein kinase Krs-2

## 产品验证



Western blot analysis of MST1/STK4 expressed in Mouse brain,Rat liver,293T using MST1/STK4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)