

# PQBP1 Rabbit pAb

货号: **AYP10431**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	6kDa/7kDa/8kDa/14kDa/17kDa/18kDa/26kDa/30kDa
实测分子量	30kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	22Rv1,SW480,HL-60,Mouse liver,Mouse brain,Rat brain,Rat ovary,Rat liver,SKOV3
细胞定位	Nucleus,Nucleus speckle
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-265 of human PQBP1 (NP_005701.1).
序列	MPLPVALQTRLAKRGILKHLEPEPEEEIIAEDYDDDPVDYEATRLEGLPPSWYKVFDPSCGLPYYWNADTDLVSWLSPHDP NSVVTKSAKKLRSSNADAEKLRSHDKSDRGHDKSDRSHEKLRGHDKSDRGHDKSDRDRERERERERERE RDRDRGYDKADREEGKERRHHRREELAPYPKSKKAVSRKDEELDPMPSSYSYDAPRGTWSTGLPKRNEAKTGADTTAA GPLFQRPYPSPGAVLRANAASRTKQQD

## 靶点信息

研究背景	This gene encodes a nuclear polyglutamine-binding protein that is involved with transcription activation. The encoded protein contains a WW domain. Mutations in this gene have been found in patients with Ren penning syndrome 1 and other syndromes with X-linked mental retardation. Multiple alternatively spliced transcript variants that encode different protein isoforms have been described for this gene.
------	-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------

基因ID	10084
基因名	PQBP1
Swiss	O60828
别名	PQBP1;MRX2;MRX55;MRXS3;MRXS8;NPW38;RENS1;SHS

## 产品验证



Western blot analysis of PQBP1 expressed in Mouse liver,Rat brain,22Rv1 using PQBP1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)