

— ABLYBIO, Help Your Research



KIAA1524 Rabbit pAb

货号: **AYP10634**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens,Mus musculus,Rattus norvegicus Co-IP: Homo sapiens,Mus musculus,Rattus norvegicus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	84kDa/102kDa
实测分子量	102kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HepG2,U-87MG,MCF7,HeLa,A375,Rat testis
细胞定位	Cytoplasm,Membrane,Single-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 681-905 of human KIAA 1524 (NP_065941.2).
------	---

靶点信息

研究背景	Oncoprotein that inhibits PP2A and stabilizes MYC in human malignancies. Promotes anchorage-independent cell growth and tumor formation.
基因ID	57650
基因名	KIAA1524
Swiss	Q8TCG1 (https://www.uniprot.org/uniprotkb/Q8TCG1/entry)
别名	KIAA1524,CIP2A,p90,KIAA1524 Rabbit pAb,Cancerous inhibitor of PP2A,p90 autoantigen

产品验证



Western blot analysis of KIAA1524 expressed in Rat testis, HepG2, U-87MG using KIAA1524 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)