

TGFB3 Rabbit pAb

货号: **AYP10690**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	35kDa/47kDa
实测分子量	60kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474,K-937,Mouse lung,Mouse skeletal muscle,Mouse heart,Rat lung,Rat heart
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 25-300 of human TGFB3 (NP_003230.1).
序列	STCTTLDFGHIKKRVEAIRGQILSKLRLTSPPEPTVMTHVPYQVLALYNSTRELLEEMHGEREEGCTQENTESEYYAKEIHKF DMIQGLAEHNELAVCPKGITSKVFRFNVSVEKNRTNLFRAEFRVLRVNPSSKRNEQRIELFQILRPDEHIAKQRYIGGKNL PTRGTAEWLSFDVTDTVREWLLRRESNLGLEISIHCPCHTFQPNGDILENIHEVMEIKFKGVDNEDDHGRGDLGRLKKQK DHHNPHLILMMIPPHRLDNPQGQGGQRKKR

靶点信息

研究背景	This gene encodes a member of the TGF-beta family of proteins. The encoded protein is secreted and is involved in embryogenesis and cell differentiation. Defects in this gene are a cause of familial arrhythmogenic right ventricular dysplasia 1.
------	--

基因ID	7043
基因名	TGFB3
Swiss	P10600
别名	TGFB3;ARVD;ARVD1;LDS5;RNHF;TGF-beta3

产品验证



Western blot analysis of TGFB3 expressed in Mouse lung,Rat lung,BT-474 using TGFB3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn