

# GPR155 Rabbit pAb

货号: **AYP10801**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	96kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 721-870 of human GPR155 (NP_001028217.1).
序列	PFKRRLEFLWNNKDTAENRDSPVSEEIKMTCQQFIHYHRDLCIRNIVKERRCGAKTSAGTFCGCDLVSWLIEVGLASDRGE AVIYGDRLVQGGVIQHITNEYEFRDEYLFYRFLQKSPEQSPAINANTLQQUERYKEIEHSSPPSHSPKT

## 靶点信息

研究背景	Involved in cognition. Located in extracellular exosome.
基因ID	151556
基因名	GPR155

Swiss	Q7Z3F1
别名	GPR155;DEP.7;DEPDC3;PGR22

## 产品验证



Western blot analysis of GPR155 expressed in Mouse brain,Rat brain using GPR155 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)