

— ABLYBIO, Help Your Research



INTS7 Rabbit pAb

货号: **AYP10840**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	101kDa/104kDa/105kDa/106kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,A-431,U-87MG,HeLa,Mouse testis
细胞定位	Chromosome,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 702-948 of human INTS 7 (NP_001186740.1).
------	---

靶点信息

研究背景	This gene encodes a subunit of the integrator complex. The integrator complex associates with the C-terminal domain of RNA polymerase II and mediates 3'-end processing of the small nuclear RNAs U1 and U2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	25896
基因名	INTS7
Swiss	Q9NVH2 (https://www.uniprot.org/uniprotkb/Q9NVH2/entry)
别名	INTS7,C1orf73,INT7,INTS7 Rabbit pAb

产品验证



Western blot analysis of INTS7 expressed in Mouse testis,Raji,A-431 using INTS7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)