

INTS7 Rabbit pAb

货号: **AYP10840**

产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	101kDa/104kDa/105kDa/106kDa
实测分子量	107kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,A-431,U-87MG,HeLa,Mouse testis
细胞定位	Chromosome,Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 702-948 of human INTS 7 (NP_001186740.1).
序列	NVELQQQSCLLISHAIEALILDPEASFQEYGSTGTAHADSEYERRMMSVYNHVLEEHTACLCNAIALLKVPLSFQRYFFQK LQSTSIKLALSPSPRNPAEPIAVQNNQLALKVEGVVQHGSKPGLFRKIQSVCLNVSSTLQSKSGQDYKIPIDNMTNEMEQ RVEPHNDYFSTQFLNFAILGTHNITVESSVKDANGIVWKTGPRTTIFVKSLEDPYSQQIRLQQQQAQQPLQQQQQRNAY TRF

靶点信息

研究背景	This gene encodes a subunit of the integrator complex. The integrator complex associates with the C-terminal domain of RNA polymerase II and mediates 3'-end processing of the small nuclear RNAs U1 and U2. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene.
基因ID	25896

基因名	INTS7
Swiss	Q9NVH2
别名	INTS7;C1orf73;INT7

产品验证



Western blot analysis of INTS7 expressed in Mouse testis,Raji,A-431 using INTS7 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn