

# NLGN1 Rabbit pAb

货号: **AYP10848**

## 产品信息

反应	Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	92kDa/93kDa
实测分子量	105kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal, 50% glycerol, pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	Cell junction, Cell membrane, Single-pass type I membrane protein, postsynaptic cell membrane, postsynaptic density, synapse
纯化	Affinity purification

## 抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 200-300 of human NLGN1 (NP_055747.1).
序列	YGNVIVITVNYRLGVLGFLSTGDQAAKGNYGLLDLIQALRWTSNIGFFGGDPLRITVFGSGAGGSCVNLLTSHYSEGNRWSNSTKGLFQRAIAQSGTAL

## 靶点信息

研究背景	This gene encodes a member of a family of neuronal cell surface proteins. Members of this family may act as splice site-specific ligands for beta-neurexins and may be involved in the formation and remodeling of central nervous system synapses.
基因ID	22871

基因名	NLGN1
Swiss	Q8N2Q7
别名	NLGN1;NL1

## 产品验证



Western blot analysis of NLGN1 expressed in Mouse kidney using NLGN1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)