

# ZNF70 Rabbit pAb

货号: **AYP10863**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	50kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse eye,Rat eye
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZNF70 (NP_068735.1).
序列	MEVPPATKFGETFAFENRLESQQGLFPGEDLGDPLFLQERGLEQMAVIYKEIPLGEQDEENDDYEGNFSLCSSPVQHQSIIP GTRPQDDELFGQTFLQKSDLSCQIIHSEEPSPCDCAETDRGDSGPNAPHRTPQPAKPY

## 靶点信息

研究背景	May be involved in transcriptional regulation.
基因ID	7621
基因名	ZNF70

Swiss	Q9UC06
别名	ZNF70;Cos17

## 产品验证



Western blot analysis of ZNF70 expressed in Jurkat, Mouse eye, Rat eye using ZNF70 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)