

— ABLYBIO, Help Your Research



ZNF70 Rabbit pAb

货号: **AYP10863**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	50kDa
实测分子量	51kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Jurkat,Mouse eye,Rat eye
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-140 of human ZNF70 (NP_068735.1).
------	---

靶点信息

研究背景	May be involved in transcriptional regulation.
------	--

基因ID	7621
基因名	ZNF70
Swiss	Q9UC06 (https://www.uniprot.org/uniprotkb/Q9UC06/entry)
别名	ZNF70,Cos17,ZNF70 Rabbit pAb,Zinc finger protein N27C7-1

产品验证



Western blot analysis of ZNF70 expressed in Jurkat, Mouse eye, Rat eye using ZNF70 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)