

# MKRN3 Rabbit pAb

货号: **AYP10864**

## 产品信息

反应	Human,Mouse
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	55kDa
实测分子量	56kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse kidney
细胞定位	
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-260 of human MKRN3 (NP_005655.1).
序列	MEEPAAPSEAHEAAGAQAQAGAEAAREGVSGPDLVCEPSGESAAPDSALPHAARGWAPFPVAPVPAHLRRGGLRPAPASG GGAWPSPLPSRSSGIWTKQIICRYIHHGQCKEGENCRYSHDLSGRKMATEGGVSPPGASAGGGPSTAAHIEPPTQEVAEA PPAASSLSLPVIGSAAERGFEEAERDNADRGAAAGGAGVESWADAIEFVPGQPYPYRGRWVASAPEAPLQSSETERKQMAVG SGLRFCYYASRGVCFRGESCMY

## 靶点信息

研究背景	The protein encoded by this gene contains a RING (C3HC4) zinc finger motif and several C3H zinc finger motifs. This gene is intronless and imprinted, with expression only from the paternal allele. Disruption of the imprinting at this locus may contribute to Prader-Willi syndrome. An antisense RNA of unknown function has been found overlapping this gene.
------	---

基因ID	7681
基因名	MKRN3
Swiss	Q13064
别名	MKRN3;CPPB2;D15S9;RNF63;ZFP127;ZNF127

## 产品验证



Western blot analysis of MKRN3 expressed in Mouse kidney using MKRN3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)