

# GBX2 Rabbit pAb

货号: **AYP10893**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:1000
理论分子量	37kDa
实测分子量	34KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG,SH-SY5Y
细胞定位	Nucleus
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-250 of human GBX2 ( NP_001476.2).
序列	MSAAFPPSLMMMQRPLGSSTAFSIDSLIGSPPQSPGHFVYTGYPMPYRPVVLPPPPPPPPALPQAALQPALPPAHPHH QIPSLPTGFCSSLAQGMALTSTLMATLPGGFSASPQHQEAARFAPQPLPGGGNFDKAEALQADAEDGKGF LLAFSAAETVQASLVGAVRGQKDESKVEDDPKGKEESFSLESVDVYSSDDNLTGQAAHKEEDPGHALEETPPSSGAAG STTSTGKNRRR

## 靶点信息

研究背景	May act as a transcription factor for cell pluripotency and differentiation in the embryo.
基因ID	2637
基因名	GBX2

Swiss	P52951
别名	GBX2

## 产品验证



Western blot analysis of GBX2 expressed in U-87MG,SH-SY5Y using GBX2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5 % non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)