

# CREBBP Rabbit pAb

货号: **AYP10898**

## 产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	265kDa
实测分子量	300kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Mouse brain,Rat brain
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-445 of mouse CREBBP (Q6GQV9).
序列	MAENLLDGPPNPKRAKLSSPGFSANDNTDFGSLFDLENDLPDELIPNGELSLNLSGNLVPDAASKHKQLSELLRGGSGSSI NPGIGNVSASSPVQQGLGGQAQGPNSTNMAASLGAMGKSPLNQGDSSTPNLPKQAASTSGPTPPASQALNPQAQKQV GLVTSSPATSQTGPGICMNANFNQTHPGLLNSNSGHSMLNQAQQGQAQVMNGSLGAAGRGRGAGMPYPAPAMQGAT SSVLAETLTQVSPQMAGHAGLNTAQAGGMTKMGMTGTTSPFGQPFSQTGGQQMGATGVNPNQLASKQSMVNSLPAFPT DIKNTSVTTVPNMSQLQTSVGVPTQAIATGPTADPEKRKLIQQQLVLLLHAHKCQRREQANGEVRACSLPHCRTMKNVLN HMTHCQAGKACQVAHCASSRQIISHWKNCTRHDPCVCLPLKNASDKRNQQTIL

## 靶点信息

研究背景	This gene is ubiquitously expressed and is involved in the transcriptional coactivation of many different transcription factors. First isolated as a nuclear protein that binds to cAMP-response element binding protein (CREB), this gene is now known to play critical roles in embryonic development, growth control, and homeostasis by coupling chromatin remodeling to transcription factor recognition. The protein encoded by this gene has intrinsic histone acetyltransferase activity and also acts as a scaffold to stabilize additional protein interactions with the transcription complex. This protein acetylates both histone and non-histone proteins. This protein shares regions of very high sequence similarity with protein p300 in its bromodomain, cysteine-histidine-rich regions, and histone acetyltransferase domain. Mutations in this gene cause Rubinstein-Taybi syndrome (RTS). Chromosomal translocations involving this gene have been associated with acute myeloid leukemia. Alternative splicing results in multiple transcript variants encoding different isoforms.
基因ID	12914
基因名	
Swiss	P45481
别名	CREBBP;CBP;KAT3A;RSTS;RSTS1

## 产品验证



Western blot analysis of CREBBP expressed in Mouse brain,Rat brain using CREBBP Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)