

— ABLYBIO, Help Your Research



MS4A8B Rabbit pAb

货号: **AYP10913**

产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	26kDa
实测分子量	26kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-87MG
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human MS4A8B (NP_113645.1).
------	---

靶点信息

研究背景	This gene encodes a member of the membrane-spanning 4A gene family. Members of this protein family are characterized by common structural features and similar intron/exon splice boundaries and display unique expression patterns among hematopoietic cells and nonlymphoid tissues. The gene encoding this protein is localized to 11q12.3, among a cluster of family members.
基因ID	83661
基因名	MS4A8
Swiss	Q9BY19 (https://www.uniprot.org/uniprotkb/Q9BY19/entry)
别名	MS4A8,4SPAN4,CD20L5,MS4A4,MS4A8B,MS4A8B Rabbit pAb,Four-span transmembrane protein 4,Membrane-spanning 4-domains subfamily A member 8B

产品验证



Western blot analysis of MS4A8B expressed in U-87MG using MS4A8B Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)