

— ABLYBIO, Help Your Research



CDKN2A/p16INK4a Rabbit pAb

货号: **AYP11013**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Mus musculus , Rattus norvegicus WB: Homo sapiens , Mus musculus , Rattus norvegicus IHC: Homo sapiens , Mus musculus , Rattus norvegicus WB IF: Mus musculus ICC: Homo sapiens,Mus musculus IHC, IF, WB: Rattus norvegicus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	8kDa/11kDa/12kDa/13kDa/16kDa/17kDa
实测分子量	16KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,HeLa
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 50-156 of human CDKN2A/p16INK4a (NP_000068.1).
------	---

靶点信息

研究背景	This gene generates several transcript variants which differ in their first exons. At least three alternatively spliced variants encoding distinct proteins have been reported, two of which encode structurally related isoforms known to function as inhibitors of CDK4 kinase. The remaining transcript includes an alternate first exon located 20 Kb upstream of the remainder of the gene; this transcript contains an alternate open reading frame (ARF) that specifies a protein which is structurally unrelated to the products of the other variants. This ARF product functions as a stabilizer of the tumor suppressor protein p53 as it can interact with, and sequester, the E3 ubiquitin-protein ligase MDM2, a protein responsible for the degradation of p53. In spite of the structural and functional differences, the CDK inhibitor isoforms and the ARF product encoded by this gene, through the regulatory roles of CDK4 and p53 in cell cycle G1 progression, share a common functionality in cell cycle G1 control. This gene is frequently mutated or deleted in a wide variety of tumors, and is known to be an important tumor suppressor gene.
基因ID	1029
基因名	CDKN2A
Swiss	P42771 (https://www.uniprot.org/uniprotkb/P42771/entry), Q8N726 (https://www.uniprot.org/uniprotkb/Q8N726/entry)
别名	CDKN2A,ARF,CDK4I,CDKN2,CMM2,INK4,INK4A,MLM,MTS-1,MTS1,P14,P14ARF,P16,P16-INK4A,P16INK4,P16INK4A,P19,P19ARF,TP16,CDKN2A/p16INK4a Rabbit pAb,Cyclin-dependent kinase 4 inhibitor A,Multiple tumor suppressor 1,Alternative reading frame,Cyclin-dependent kinase inhibitor 2A

产品验证



Western blot analysis of CDKN2A/p16INK4a expressed in 293T, HeLa using CDKN2A/p16INK4a Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)