

— ABLYBIO, Help Your Research



Phospho-p70 S6 Kinase 1-T389 Rabbit pAb

货号: **AYP11049**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: MSC cells , unknow , Homo sapiens , Mus musculus , Danio rerio , Sus scrofa , M. japonicus , Rattus norvegicus , Spodoptera frugiperda , Gallus gallus IHC: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:100 - 1:500 IF/ICC: 1:50 - 1:200
理论分子量	51kDa/52kDa/56kDa/59kDa
实测分子量	70KDa/85KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	NIH/3T3
细胞定位	Cell junction,Cytoplasm,Mitochondrion outer membrane,Mitochondrion,Nucleus,synapse,synaptosome
纯化	Affinity purification


抗原信息

抗原信息	A synthetic phosphorylated peptide around T389 of human P70 S6K (NP_003152.1).
------	--

靶点信息

研究背景	This gene encodes a member of the ribosomal S6 kinase family of serine/threonine kinases. The encoded protein responds to mTOR (mammalian target of rapamycin) signaling to promote protein synthesis, cell growth, and cell proliferation. Activity of this gene has been associated with human cancer. Alternatively spliced transcript variants have been observed. The use of alternative translation start sites results in isoforms with longer or shorter N-termini which may differ in their subcellular localizations. There are two pseudogenes for this gene on chromosome 17.
基因ID	6198
基因名	RPS6KB1
Swiss	P23443 (https://www.uniprot.org/uniprotkb/P23443/entry)
别名	PS6K,S6K,S6K-beta-1,S6K1,STK14A,p70 S6KA,p70(S6K)-alpha,p70-S6K,p70-alpha,P70 S6K,RPS6KB1,p70S6KA,Phospho-p70 S6 Kinase 1-T389 Rabbit pAb,70 kDa ribosomal protein S6 kinase 1,Ribosomal protein S6 kinase 1,Serine/threonine-protein kinase 14A,p70 ribosomal S6 kinase alpha

产品验证



Western blot analysis of Phospho-p70 S6 Kinase 1-T389 expressed in NIH/3T3 using Phospho-p70 S6 Kinase 1-T389 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)