

# CD133 Rabbit pAb

货号: **AYP11219**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>IF:</b> Homo sapiens , Mus musculus <b>IHC:</b> Homo sapiens , Mus musculus <b>WB:</b> Homo sapiens
应用	<a href="#">WB</a> <a href="#">IHC</a> <a href="#">IF/ICC</a> <a href="#">IP</a>
推荐浓度	<b>WB:</b> 1:500 - 1:1000 <b>IHC:</b> 1:50 - 1:200 <b>IF/ICC:</b> 1:50 - 1:200 <b>IP:</b> 1:50 - 1:200
理论分子量	92kDa/93kDa/94kDa/96kDa/97kDa
实测分子量	120KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HCT116,Mouse brain,Rat lung
细胞定位	Apical cell membrane,Cell projection,Endoplasmic reticulum,Endoplasmic reticulum-Golgi intermediate compartment,Multi-pass membrane protein,cilium,microvillus membrane,photoreceptor outer segment
纯化	Affinity purification


## 抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 750-850 of human CD133 (NP_006008.1).
序列	FEHYLQWIEFSISEKVASCKPVATALDTAVDVFVLCVSIIDPLNLFWFGIGKATVFLPALIFAVKLAKEYRMRMDSSEVYDDVETIPMKNMENGNGYHKDH

## 靶点信息

研究背景	This gene encodes a pentaspan transmembrane glycoprotein. The protein localizes to membrane protrusions and is often expressed on adult stem cells, where it is thought to function in maintaining stem cell properties by suppressing differentiation. Mutations in this gene have been shown to result in retinitis pigmentosa and Stargardt disease. Expression of this gene is also associated with several types of cancer. This gene is expressed from at least five alternative promoters that are expressed in a tissue-dependent manner. Multiple transcript variants encoding different isoforms have been found for this gene.
基因ID	8842
基因名	PROM1
Swiss	O43490
别名	PROM1;AC133;CD133;CORD12;MCDR2;MSTP061;PROML1;RP41;STGD4

## 产品验证



Western blot analysis of CD133 expressed in HCT116, Mouse brain, Rat lung using CD133 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)