

Nestin Rabbit pAb

货号: **AYP11227**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mouse Monocyte Macrophage Leukemia Cells , Homo sapiens , Mus musculus IHC: Mus musculus , Rattus norvegicus IF: Mus musculus , Rattus norvegicus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	177kDa
实测分子量	150KDa/177KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Neuro-2a
细胞定位	cytoplasm,intermediate filament cytoskeleton
纯化	Affinity purification

抗原信息


抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1522-1621 of human NES (NP_006608.1).
序列	GMEDAGPGADIIGVNGQGPNLEGKSQLVNGGVMNGLEQSEEVGQGMPLVSEGDRGSPFQEEEGSALKTSWAGAPVH LGQQQLKFTQREGDRESWSSGED

靶点信息

研究背景	This gene encodes a member of the intermediate filament protein family and is expressed primarily in nerve cells.
------	---

基因ID	10763
基因名	NES
Swiss	P48681
别名	NES;Nbla00170;nestin;Nestin

产品验证



Western blot analysis of Nestin expressed in Neuro-2a using Nestin Rabbit pAb at 1:1000. Secondary anti body: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% no n-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn