

— ABLYBIO, Help Your Research



FTO Rabbit pAb

货号: **AYP11277**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Homo sapiens , Rattus norvegicus IHC: Mus musculus
应用	WB IP
推荐浓度	WB: 1:500 - 1:1000 IP: 1:500 - 1:1000
理论分子量	6kDa/12kDa/14kDa/58kDa
实测分子量	60KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,SH-SY5Y
细胞定位	Nucleus,Nucleus speckle
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human FTO (N P_001073901.1).
------	---

靶点信息

研究背景	This gene is a nuclear protein of the AlkB related non-haem iron and 2-oxoglutarate-dependent oxygenase superfamily but the exact physiological function of this gene is not known. Other non-heme iron enzymes function to reverse alkylated DNA and RNA damage by oxidative demethylation. Studies in mice and humans indicate a role in nervous and cardiovascular systems and a strong association with body mass index, obesity risk, and type 2 diabetes.
基因ID	79068
基因名	FTO
Swiss	Q9C0B1 (https://www.uniprot.org/uniprotkb/Q9C0B1/entry)
别名	FTO,ALKBH9,BMIQ14,GDFD,FTO Rabbit pAb,Fat mass and obesity-associated protein,U6 small nuclear RNA (2'-O-methyladenosine-N(6)-)-demethylase FTO,U6 small nuclear RNA N(6)-methyladenosine-demethylase FTO,mRNA (2'-O-methyladenosine-N(6)-)-demethylase FTO

产品验证



Western blot analysis of FTO expressed in HeLa,293T,SH-SY5Y using FTO Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)