

— ABLYBIO, Help Your Research



MT-CYB Rabbit pAb

货号: **AYP11301**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mouse heart , Homo sapiens , Danio rerio , Mus musculus , Sus scrofa
应用	WB
推荐浓度	WB: 1:200 - 1:1000
理论分子量	43kDa
实测分子量	43kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Raji,Rat testis
细胞定位	Mitochondrion inner membrane,Multi-pass membrane protein
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1-100 of mouse MT-CYB (NP_904340.1).
------	---

靶点信息

研究背景	Component of the ubiquinol-cytochrome c reductase complex (complex III or cytochrome b-c1 complex) that is part of the mitochondrial respiratory chain. The b-c1 complex mediates electron transfer from ubiquinol to cytochrome c. Contributes to the generation of a proton gradient across the mitochondrial membrane that is then used for ATP synthesis.
基因ID	17711
基因名	Mt-Cyb
Swiss	P00158
别名	MT-CYB,MTCYB,CYTB,MT-CYB Rabbit pAb,Complex III subunit 3,Complex III subunit III,Cytochrome b-c1 complex subunit 3,Ubiquinol-cytochrome-c reductase complex cytochrome b subunit,Cob,mt-Cytb

产品验证



Western blot analysis of MT-CYB expressed in Raji, Rat testis using MT-CYB Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)