

SIRT3 Rabbit pAb

货号: **AYP11302**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: 769-P PLVX;769-P SIRT3 OE , Mus musculus , Homo sapiens , Rattus norvegicus IF: Homo sapiens , Mus musculus IHC: Homo sapiens , Mus musculus
应用	WB IHC IP IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IP: 1:50 - 1:500 IF/ICC: 1:50 - 1:200
理论分子量	28kDa/43kDa
实测分子量	28KDa/44KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,293T,HepG2,Mouse kidney
细胞定位	Mitochondrion matrix
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 300-399 of human SIRT3 (NP_03637 1.1).
序列	QRFLLVVDFPMADLLLLILGTSLEVEPFASLTEAVRSSVPRLLINRDLVGPLAWHPRSRDVAQLGDVVHGVESLVELLGWT EEMRDLVQRETGKLDGPKD

靶点信息

研究背景	This gene encodes a member of the sirtuin family of proteins, homologs to the yeast Sir2 protein. Members of the sirtuin family are characterized by a sirtuin core domain and grouped into four classes. The functions of human sirtuins have not yet been determined; however, yeast sirtuin proteins are known to regulate epigenetic gene silencing and suppress recombination of rDNA. Studies suggest that the human sirtuins may function as intracellular regulatory proteins with mono-ADP-ribosyltransferase activity. The protein encoded by this gene is included in class I of the sirtuin family. Two alternatively spliced transcript variants that encode different proteins have been described for this gene.
基因ID	23410
基因名	SIRT3
Swiss	Q9NTG7
别名	SIRT3;SIR2L3;sirtuin 3

产品验证



Western blot analysis of SIRT3 expressed in Mouse kidney, HeLa, 293T using SIRT3 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn