

— ABLYBIO, Help Your Research



ITGAX Rabbit pAb

货号: **AYP11383**

产品信息

| | |
|-------|--|
| 反应 | Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | IHC: Mus musculus , Rattus norvegicus , Sus scrofa IF: Homo sapiens,Mus musculus , Rattus norvegicus WB: Mus musculus |
| 应用 | IF/ICC |
| 推荐浓度 | IF/ICC: 1:50 - 1:200 |
| 理论分子量 | 127kDa |
| 实测分子量 | 127kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | |
| 细胞定位 | Membrane,Single-pass type I membrane protein |
| 纯化 | Affinity purification |

抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 146-342 of ITGAX (NP_00878.2). |
|------|--|

靶点信息

| | |
|-------|---|
| 研究背景 | This gene encodes the integrin alpha X chain protein. Integrins are heterodimeric integral membrane proteins composed of an alpha chain and a beta chain. This protein combines with the beta 2 chain (ITGB2) to form a leukocyte-specific integrin referred to as inactivated-C3b (iC3b) receptor 4 (CR4). The alpha X beta 2 complex seems to overlap the properties of the alpha M beta 2 integrin in the adherence of neutrophils and monocytes to stimulated endothelium cells, and in the phagocytosis of complement coated particles. Two transcript variants encoding different isoforms have been found for this gene. |
| 基因ID | 3687 |
| 基因名 | ITGAX |
| Swiss | P20702 (https://www.uniprot.org/uniprotkb/P20702/entry) |
| 别名 | ITGAX,CD11C,SLEB6,ITGAX Rabbit pAb,CD11 antigen-like family member C,Leu M5,Leukocyte adhesion glycoprotein p150,95 alpha chain,Leukocyte adhesion receptor p150,95 |

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)