



PER2 Rabbit pAb

货号: **AYP11585**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus , Rattus norvegicus IHC: Homo sapiens,Mus musculus , Rattus norvegicus
应用	WB IHC ChIP
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 ChIP: 1:50 - 1:200
理论分子量	45kDa/136kDa
实测分子量	136KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse liver
细胞定位	Cytoplasm,Nucleus,nucleolus,perinuclear region
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 390-565 of human PER2 (NP_073728.1).
------	--

靶点信息

研究背景	This gene is a member of the Period family of genes and is expressed in a circadian pattern in the suprachiasmatic nucleus, the primary circadian pacemaker in the mammalian brain. Genes in this family encode components of the circadian rhythms of locomotor activity, metabolism, and behavior. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene may increase the risk of getting certain cancers and have been linked to sleep disorders.
基因ID	8864
基因名	PER2
Swiss	O15055 (https://www.uniprot.org/uniprotkb/O15055/entry)
别名	PER2,FASPS,FASPS1,PER2 Rabbit pAb,Circadian clock protein PERIOD 2,KIAA0347

产品验证



实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)