

— ABLYBIO, Help Your Research



Pea3 / ETV4 Rabbit pAb

货号: **AYP11689**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens WB: Mus musculus , Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	23kDa/49kDa/53kDa
实测分子量	54KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,SW480,Mouse stomach,Mouse heart,Rat heart
细胞定位	Nucleus
纯化	Affinity purification


抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 385-484 of human ETV4 (NP_001977.1).
------	---

靶点信息

研究背景	Transcriptional activator. May play a role in keratinocyte differentiation.
基因ID	2118
基因名	ETV4
Swiss	P43268 (https://www.uniprot.org/uniprotkb/P43268/entry)
别名	ETV4,E1A-F,E1AF,PEA3,PEAS3,ETS variant 4,Pea3/ETV4,Pea3 / ETV4 Rabbit pAb,Adenovirus E1A enhancer-binding protein,Polyomavirus enhancer activator 3 homolog

产品验证



Western blot analysis of Pea3 / ETV4 expressed in Mouse stomach,Rat heart,HeLa using Pea3 / ETV4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)