

— ABLYBIO, Help Your Research



ABCB11 Rabbit pAb

货号: **AYP11862**

产品信息

反应	Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IF: Homo sapiens WB: Homo sapiens,Mus musculus
应用	IF/ICC
推荐浓度	IF/ICC: 1:50 - 1:200
理论分子量	146kDa
实测分子量	146kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	Rat liver
细胞定位	Membrane,Multi-pass membrane protein
纯化	Affinity purification

抗原信息

抗原信息	A synthetic peptide corresponding to a sequence within amino acids 1050-1150 of human ABCB11 (NP_03733.2).
------	--

靶点信息

研究背景	The membrane-associated protein encoded by this gene is a member of the superfamily of ATP-binding cassette (ABC) transporters. ABC proteins transport various molecules across extra- and intra-cellular membranes. ABC genes are divided into seven distinct subfamilies (ABC1, MDR/TAP, MRP, ALD, OABP, GCN20, White). This protein is a member of the MDR/TAP subfamily. Members of the MDR/TAP subfamily are involved in multidrug resistance. The protein encoded by this gene is the major canalicular bile salt export pump in man. Mutations in this gene cause a form of progressive familial intrahepatic cholestases which are a group of inherited disorders with severe cholestatic liver disease from early infancy.
基因ID	8647
基因名	ABCB11
Swiss	O95342 (https://www.uniprot.org/uniprotkb/O95342/entry)
别名	ABCB11,ABC16,BRIC2,BSEP,PFIC-2,PFIC2,PGY4,SPGP,ABCB11 Rabbit pAb,ATP-binding cassette sub-family B member 11

产品验证

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)