

# PLGF Rabbit pAb

货号: **AYP11985**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Rattus norvegicus , Homo sapiens
应用	<a href="#">WB</a>
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	16kDa/19kDa/24kDa/27kDa
实测分子量	25kDa/50kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	293T,Mouse brain,Mouse liver,HT-29,MCF7,Rat brain
细胞定位	Secreted
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 19-170 of human PLGF ( NP_002623.2).
序列	LPAVPPQWALSAGNGSSEVEVVPFQEVWGRSYCRALERLVDVVSEYPSVEHEMFSPSCVSLLRCTGCCGDENLHCVPV ETANVTMQLLKIRSGDRPSYVELTFSQHVRCECRPLREKMKPERRRPKGRGKRRREKQRPTDCHLCGDAVPRR

## 靶点信息

研究背景	This gene encodes a growth factor found in placenta which is homologous to vascular endothelial growth factor. Alternatively spliced transcripts encoding different isoforms have been found for this gene.
基因ID	5228
基因名	PGF

<b>Swiss</b>	P49763
别名	PGF;D12S1900;PGFL;PIGF;PLGF;PIGF-2;SHGC-10760

## 产品验证



Western blot analysis of PLGF expressed in Mouse brain,Rat brain,293T using PLGF Rabbit pAb at 1:100. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn)