

— ABLYBIO, Help Your Research



# SERPINB2 Rabbit pAb

货号: **AYP11998**

## 产品信息

反应	Human
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Homo sapiens , Rattus norvegicus
应用	WB
推荐浓度	<b>WB:</b> 1:200 - 1:2000
理论分子量	46kDa
实测分子量	47kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	U-251MG,Jurkat,LO2
细胞定位	Cytoplasm,Secreted,extracellular space
纯化	Affinity purification

## 抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-415 of human SERPIN B2 (NP_002566.1).
------	---

## 靶点信息

研究背景	Inhibits urokinase-type plasminogen activator. The monocyte derived PAI-2 is distinct from the endothelial cell-derived PAI-1.
------	--

基因ID	5055
基因名	SERPINB2
Swiss	P05120 ( <a href="https://www.uniprot.org/uniprotkb/P05120/entry">https://www.uniprot.org/uniprotkb/P05120/entry</a> )
别名	SERPINB2,HsT1201,PAI,PAI-2,PAI2,PLANH2,SERPINB2 Rabbit pAb,Monocyte Arg-serpin,Placental plasminogen activator inhibitor,Serpin B2,Urokinase inhibitor

## 产品验证



Western blot analysis of SERPINB2 expressed in U-251MG, Jurkat, LO2 using SERPINB2 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/www.ablybio.cn>)