


研究背景	The protein encoded by this gene is a 47 kDa cytosolic subunit of neutrophil NADPH oxidase. This oxidase is a multicomponent enzyme that is activated to produce superoxide anion. Mutations in this gene have been associated with chronic granulomatous disease.
基因ID	653361
基因名	NCF1
Swiss	P14598
别名	NCF1;NCF1A;NOXO2;SH3PXD1A;p47phox

产品验证



Western blot analysis of p47phox/NCF1 expressed in Raji, Daudi, Jurkat (Negative control) using p47phox/NCF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn