

— ABLYBIO, Help Your Research



IFNAR1 Rabbit pAb

货号: **AYP12118**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	IHC: Homo sapiens WB IF: Mus musculus
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	49kDa/56kDa/63kDa
实测分子量	64KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	SP2/0,HeLa,A-549
细胞定位	Membrane,Single-pass type I membrane protein
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 470-557 of human IFNAR1 (NP_000620.2).
------	--

靶点信息

研究背景	The protein encoded by this gene is a type I membrane protein that forms one of the two chains of a receptor for interferons alpha and beta. Binding and activation of the receptor stimulates Janus protein kinases, which in turn phosphorylate several proteins, including STAT1 and STAT2. The encoded protein also functions as an antiviral factor.
基因ID	3454
基因名	IFNAR1
Swiss	P17181 (https://www.uniprot.org/uniprotkb/P17181/entry)
别名	IFNAR1,AVP,IFN-alpha-REC,IFNAR,IFNBR,IFRC,IFNAR1 Rabbit pAb,Cytokine receptor class-II member 1,Cytokine receptor family 2 member 1,Type I interferon receptor 1

产品验证



Western blot analysis of IFNAR1 expressed in SP2/0, HeLa, A-549 using IFNAR1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)