

MZF1 Rabbit pAb

货号: **AYP12235**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | WB: Mus musculus , Rattus norvegicus IHC: Rattus norvegicus |
| 应用 | WB |
| 推荐浓度 | WB: 1:500 - 1:2000 |
| 理论分子量 | 31kDa/54kDa/82kDa |
| 实测分子量 | 110kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa,A-431,U-251MG,293T,SGC-7901,Mouse brain,Mouse lung,Rat brain |
| 细胞定位 | Nucleus |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 100-350 of human MZF 1 (NP_003413.2). |
| 序列 | IQARVQGQRPGSPEEAAALVDGLRREPGGPRRWTVQVQGQEVLEKMEPSSFQPLPETEPPTPEPGPKTPPRTMQESPLGLQVKEESEVTEDESDFLESGPLAATQESVPTLLPEEAQRCGTVLDDQIFPHSKTGPEGPSWREHPRALWHEEAGGIFSPGFA LQLGSISAGPGSVSPHLHVPWDLGMAGLSGQIQSPSREGGF AHALLLPSDLRSEQDPTDEDPCRGVGPALITTRWRSRGRSRGRPSTGG |

靶点信息

| | |
|------|--|
| 研究背景 | Binds to target promoter DNA and functions as transcription regulator. Regulates transcription from the P ADI1 and CDH2 promoter. May be one regulator of transcriptional events during hemopoietic development. |
|------|--|

| | |
|-------|-------------------------------------|
| 基因ID | 7593 |
| 基因名 | MZF1 |
| Swiss | P28698 |
| 别名 | MZF1;MZF-1;MZF1B;ZFP98;ZNF42;ZSCAN6 |

产品验证



Western blot analysis of MZF1 expressed in Mouse brain,Rat brain,HeLa using MZF1 Rabbit pAb at 1:10 00. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Bl ocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn