

— ABLYBIO, Help Your Research



# MZF1 Rabbit pAb

货号: **AYP12235**

## 产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	<b>WB:</b> Mus musculus , Rattus norvegicus <b>IHC:</b> Rattus norvegicus
应用	WB
推荐浓度	<b>WB:</b> 1:500 - 1:2000
理论分子量	31kDa/54kDa/82kDa
实测分子量	110kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,A-431,U-251MG,293T,SGC-7901,Mouse brain,Mouse lung,Rat brain
细胞定位	Nucleus
纯化	Affinity purification


## 抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 100-350 of human MZF1 (NP_003413.2).
------	--

## 靶点信息

研究背景	Binds to target promoter DNA and functions as transcription regulator. Regulates transcription from the P ADI1 and CDH2 promoter. May be one regulator of transcriptional events during hemopoietic development.
基因ID	7593
基因名	MZF1
Swiss	P28698 ( <a href="https://www.uniprot.org/uniprotkb/P28698/entry">https://www.uniprot.org/uniprotkb/P28698/entry</a> )
别名	MZF1,MZF-1,MZF1B,ZFP98,ZNF42,ZSCAN6,MZF1 Rabbit pAb,Zinc finger and SCAN domain-containing protein 6,Zinc finger protein 42,MZF

## 产品验证



Western blot analysis of MZF1 expressed in Mouse brain, Rat brain, HeLa using MZF1 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

## 实验步骤

访问官网浏览详情: [www.ablybio.cn](http://www.ablybio.cn) (<https://www.ablybio.cn/>[www.ablybio.cn](http://www.ablybio.cn))