

ISG15 Rabbit pAb

货号: **AYP12262**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Dragon Li cats , Bos taurus
应用	WB IHC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200
理论分子量	17kDa
实测分子量	15KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	RD,Jurkat
细胞定位	Cytoplasm,Secreted
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-165 of human ISG15 (NP_005092.1).
序列	MGWDLTVKMLAGNEFQVSLSSSMSVSELKAQITQKIGVHAFQQLAVHPSGVALQDRVPLASQGLGPGSTVLLVVDKC DEPLSILVRNNKGRSSTYEVRLTQTVAHLKQVSGLEGVQDDLFWLTFEGKPLEDQLPLGEYGLKPLSTVFMNLRRLRGGG TEPGGRS

靶点信息

研究背景	The protein encoded by this gene is a ubiquitin-like protein that is conjugated to intracellular target proteins upon activation by interferon-alpha and interferon-beta. Several functions have been ascribed to the encoded protein, including chemotactic activity towards neutrophils, direction of ligated target proteins to intermediate filaments, cell-to-cell signaling, and antiviral activity during viral infections. While conjugates of this protein have been found to be noncovalently attached to intermediate filaments, this protein is sometimes secreted.
基因ID	9636
基因名	ISG15
Swiss	P05161
别名	ISG15;G1P2;IFI15;IMD38;IP17;UCRP;hUCRP

产品验证



Western blot analysis of ISG15 expressed in RD, Jurkat using ISG15 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn