

NUMB Rabbit pAb

货号: **AYP12271**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 IF/ICC: 1:50 - 1:200
理论分子量	37kDa/48-70kDa
实测分子量	71KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,NIH/3T3,SH-SY5Y,Mouse brain,Rat testis,Rat kidney
细胞定位	Membrane,Peripheral membrane protein
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 521-640 of human NUMB (NP_003735.3).
序列	SQMVANVFGTAGHPQAAHPHQSPSLVRQQTFFHYEASSATTSPFFKPPAQLNGSAAFNGVDDGRLASADRHTEVPTG TCPVDPFEAQWAALENKSKQRTNPSPTNPFSSDLQKTFEIEL

靶点信息

研究背景	The protein encoded by this gene plays a role in the determination of cell fates during development. The encoded protein, whose degradation is induced in a proteasome-dependent manner by MDM2, is a membrane-bound protein that has been shown to associate with EPS15, LNX1, and NOTCH1. Alternative splicing results in multiple transcript variants.
------	---

基因ID	8650
基因名	NUMB
Swiss	P49757
别名	NUMB;C14orf41;S171;c14_5527

产品验证



Western blot analysis of NUMB expressed in Mouse brain,Rat testis,HeLa using NUMB Rabbit pAb at 1:1 000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. B locking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn