

DUSP10 Rabbit pAb

货号: **AYP12273**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens
应用	WB IHC
推荐浓度	WB: 1:500 - 1:2000 IHC: 1:50 - 1:200
理论分子量	16kDa/52kDa
实测分子量	53kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	BT-474
细胞定位	Cytoplasm,Nucleus
纯化	Affinity purification


抗原信息

抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human DUSP10 (NP_009138.1).
序列	MPPSPLDDRIVVALSRPVRPQDLNLCLDSSYLGSANPGSNHPPVIATTVVSLKAANLTYMPSSSGSARSLNCGCSSASC CTVATYDKDNQAQTQAIAAGTTTTAIGTSTTCPANQMVNNNENTGSLSPSSGSPVSGTPKQLASIKIIPNDLAKKMTK CSKSHLPSQGPVIIDCRPFMEYKSHIQGAVHINCADKISRRRLQGGKITVLDLISCREGKDSFKRIFSKEIIVYDENTNEPSR VMPSQPLHIVLES�KREGKEPLVLKGGLSFFKQNHENLCDNSLQLQECREVGGA

靶点信息

研究背景	Dual specificity protein phosphatases inactivate their target kinases by dephosphorylating both the phosphoserine/threonine and phosphotyrosine residues. They negatively regulate members of the MAP kinase superfamily, which is associated with cellular proliferation and differentiation. Different members of this family of dual specificity phosphatases show distinct substrate specificities for MAP kinases, different tissue distribution and subcellular localization, and different modes of expression induction by extracellular stimuli. This gene product binds to and inactivates p38 and SAPK/JNK. Alternative splicing results in multiple transcript variants.
基因ID	11221
基因名	DUSP10
Swiss	Q9Y6W6
别名	DUSP10;MKP-5;MKP5

产品验证



Western blot analysis of DUSP10 expressed in BT-474 using DUSP10 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn