

USP26 Rabbit pAb

货号: **AYP12287**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Homo sapiens , Mus musculus
应用	WB
推荐浓度	WB: 1:500 - 1:2000
理论分子量	104kDa
实测分子量	104kDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HeLa,U-937,Mouse testis,Mouse liver,Mouse heart,Rat liver
细胞定位	Nucleus
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 664-913 of human USP26 (NP_114113.1).
序列	TSLCQFHKAGGKPASSPGTPLSKVDFQTPENPKRKKYVKTSKFVAFDRIINPTKDLYEDKNIRIPERFQKVSEQTQQCDGMRICEQAPQQALPQSFPKPGTQGHTKNLLRPTKLNLLQKSNRNSLLALGSNKNPRNKDILDKIKSKAKETKRNDKGDHTYRLISVVSHLGKTLKSGHYICDAYDFEKQIWFYDDMRVLGIQEAQMEDRRCTGYIFFYMHNEIFEEMLKREENAQLNSKEVEETLQKE

靶点信息

研究背景	This gene encodes a member of the ubiquitin-specific processing (UBP) family of proteases and is a deubiquitinating enzyme (DUB) with His and Cys domains. It is specifically expressed in testis tissue. Mutations in this gene have been associated with Sertoli cell-only syndrome and male infertility.
基因ID	83844

基因名	USP26
Swiss	Q9BXU7
别名	USP26

产品验证



Western blot analysis of USP26 expressed in Mouse testis, Rat liver, HeLa using USP26 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn