

— ABLYBIO, Help Your Research



IL22 Rabbit pAb

货号: **AYP12347**

产品信息

反应	Human,Mouse,Rat
宿主	Rabbit
克隆性	Polyclonal
预测反应	WB: Mus musculus
应用	WB IF/ICC
推荐浓度	WB: 1:500 - 1:1000 IF/ICC: 1:50 - 1:200
理论分子量	20kDa
实测分子量	16KDa
形式	Liquid
保存条件	Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.05% proclin300,50% glycerol,pH7.3.
偶联物	Unconjugated
阳性对照	HL-60
细胞定位	Secreted
纯化	Affinity purification

抗原信息


抗原信息	Recombinant fusion protein containing a sequence corresponding to amino acids 34-179 of human IL22 (NP_065386.1).
------	---

靶点信息

研究背景	Cytokine that contributes to the inflammatory response in vivo.
------	---

基因ID	50616
基因名	IL22
Swiss	Q9GZX6 (https://www.uniprot.org/uniprotkb/Q9GZX6/entry)
别名	IL22,IL-21,IL-22,IL-D110,IL-TIF,ILTIF,TIFIL-23,TIFa,zcyto18,IL22 Rabbit pAb,Cytokine Zcyto18,IL-10-related T-cell-derived-inducible factor

产品验证



Western blot analysis of IL22 expressed in HL-60 using IL22 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/www.ablybio.cn>)