

— ABLYBIO, Help Your Research



MOB4 Rabbit pAb

货号: **AYP12398**

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse,Rat |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | WB: Homo sapiens WB, IP: Homo sapiens |
| 应用 | WB IHC IF/ICC |
| 推荐浓度 | WB: 1:500 - 1:2000 IHC: 1:50 - 1:100 IF/ICC: 1:50 - 1:100 |
| 理论分子量 | 22kDa/23kDa/26kDa |
| 实测分子量 | 26kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.01% thiomersal,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | A-549,MCF7,Mouse brain,Mouse heart,Mouse kidney,Rat spleen,Rat brain |
| 细胞定位 | Cytoplasm,Golgi apparatus,Golgi stack membrane,Membrane,Peripheral membrane protein,perinuclear region |
| 纯化 | Affinity purification |


抗原信息

| | |
|------|--|
| 抗原信息 | Recombinant fusion protein containing a sequence corresponding to amino acids 1-225 of human MOB4 (NP_056202.2). |
|------|--|

靶点信息

| | |
|-------|--|
| 研究背景 | This gene was identified based on its similarity with the mouse counterpart. Studies of the mouse counterpart suggest that the expression of this gene may be regulated during oocyte maturation and preimplantation following zygotic gene activation. Alternatively spliced transcript variants encoding distinct isoforms have been observed. Naturally occurring read-through transcription occurs between this locus and the neighboring locus HSPE1. |
| 基因ID | 25843 |
| 基因名 | MOB4 |
| Swiss | Q9Y3A3 |
| 别名 | MOB4,2C4D,CGI-95,MOB1,MOB3,MOBKL3,PHOCN,PREI3,MOB4 Rabbit pAb,Class II mMOB1,Mob1 homolog 3,Mps one binder kinase activator-like 3,Preimplantation protein 3 |

产品验证



Western blot analysis of MOB4 expressed in Mouse brain,Rat spleen,A-549 using MOB4 Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn (<https://www.ablybio.cn/>www.ablybio.cn)