

FAM3A Rabbit pAb

货号: AYP12426

产品信息

| | |
|-------|--|
| 反应 | Human,Mouse |
| 宿主 | Rabbit |
| 克隆性 | Polyclonal |
| 预测反应 | WB: Mus musculus |
| 应用 | WB IHC |
| 推荐浓度 | WB: 1:500 - 1:1000 IHC: 1:50 - 1:200 |
| 理论分子量 | 21kDa/23kDa/25kDa |
| 实测分子量 | 25kDa |
| 形式 | Liquid |
| 保存条件 | Store at -20°C. Avoid freeze / thaw cycles. Buffer: PBS with 0.02% sodium azide,50% glycerol,pH7.3. |
| 偶联物 | Unconjugated |
| 阳性对照 | HeLa |
| 细胞定位 | Secreted |
| 纯化 | Affinity purification |

抗原信息


| | |
|------|------------------------------------|
| 抗原信息 | Recombinant protein of human FAM3A |
| 序列 | |

靶点信息

| | |
|------|---|
| 研究背景 | This gene encodes a cytokine-like protein. The expression of this gene may be regulated by peroxisome proliferator-activated receptor gamma, and the encoded protein may be involved in the regulation of glucose and lipid metabolism. Alternative splicing results in multiple transcript variants. |
| 基因ID | 60343 |
| 基因名 | FAM3A |

| | |
|--------------|------------------------------|
| Swiss | P98173 |
| 别名 | FAM3A;2.19;DLD;DXS560S;XAP-7 |

产品验证



Western blot analysis of FAM3A expressed in HeLa using FAM3A Rabbit pAb at 1:1000. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:5000. Lysates/proteins: 30ug per lane. Blocking buffer: 5% non-fat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 120s.

实验步骤

访问官网浏览详情: www.ablybio.cn